



BETTER MICROPLATES MEAN BETTER RESULTS

Drawing on many years of experience and in-depth knowledge, PerkinElmer offers microplates for almost every application. They are the same microplates we use to develop and validate our instruments and reagents, providing optimal results from your Alpha, TR-FRET, fluorescence, luminescence, absorbance, microfluidic, radiometric, or phenotypic imaging assays.

We also provide a range of plates and other consumables for automated liquid handling, sample prep workstations and storage. All of our microplates are manufactured using the highest quality plastics. They deliver important benefits such as low background, optimum light transmission and little or no crosstalk. We also provide custom coating, barcoding and other specialized services to meet your needs.

CulturPlates™ are sterile, Tissue Culture (TC)-treated, gamma irradiated and recommended for use with adherent or suspension cells.

The optimized well design and high-quality plastics provide low background and crosstalk. They are solid, opaque plates that are available in black or white and must be read in top-reading plate readers. They are available in 24-, 96-, 384- and 1536- well formats.

White microplates with extremely high light reflection efficiency are ideal solutions for luminescence assays, black microplates with ultra-low background fluorescence are suitable for cell-based fluorescence assays.



Exciting Offer on our entire range of culturplates™

Cat. No.	Description / Specifications	Pack size	Promo Price* (€)
6005688	culturplate-96, tc- treated / White, opaque, sterile, 96-well	160	571.20
6005689	culturplate-96, tc- treated / White, opaque, sterile, 96-well	200	547.20
6007680	culturplate-384, tc- treated / White, opaque, sterile, 384-well	50	372.00
6007689	culturplate-384, tc- treated / White, opaque, sterile, 384-well	200	1,081.20
6004680	culturplate-1536, tc- treated / White, opaque, sterile, 1536-well	50	1,147.80
6005660	culturplate-96, tc- treated / black, opaque, sterile, 96-well	50	188.40
6005669	culturplate-96, tc- treated / black, opaque, sterile, 96-well	200	547.20
6007660	culturplate-384, tc- treated / black, opaque, sterile, 384-well	50	372.00
6004660	culturplate-1536, tc- treated / black, opaque, sterile, 1536-well	50	1,146.60

Offer valid until December 31, 2022

*VAT will be added to the mentioned prices

For more information contact:

Agnes Vaino

Product manager | Nova natura

mob.: +372 53477200

email: agnes@novabio.ee



PerkinElmer, Inc.
940 Winter street
Waltham, ma 02451 usa
p: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright ©2022, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.