

Bio PP-4, platform

DESCRIPTION

Flat platform with non-slip silicone mat for Petri dishes, culture flasks, agglutination.

SPECIFICATIONS

Working area

230 × 230 mm



CAT. NUMBER

BS-010116-AK

Platform Bio PP-4

ACCESSORIES



PSU-10i
Without platform
Orbital Shaker

Shaker **PSU-10i** provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of five (5) ...

[read more](#)



CPS-20
Without platform
CO₂ Shaker

CO₂ Shaker **CPS-20** provides regulated orbital motion of the platform and is designed for use in Biosan's [S-Bt](#) CO₂ Incubator. CPS-20 is specifically designed for use in harsh ...

[read more](#)

