

SR-64, rotor

DESCRIPTION

Rotor for 8 x 8-section 0,2 ml microtube strips - for any type of strips including paired

SPECIFICATIONS

8-section strips	8
Amount of 0.2 ml tubes	64



CAT. NUMBER

BS-010201-EK

Rotor SR-64



FV-2400 Micro-Spin
Including rotors R-1.5M, R-0.5/0.2M
Mini-Centrifuge/Vortex

Mini-centrifuge/vortex Micro-Spin **FV-2400** is specially designed for genetic engineering research (for PCR-diagnostics experiments). Units can be used in microbiological, biochemical, clinical laboratories and industrial biotechnological laboratories.

Micro-Spin provides simultaneous mixing ...

[read more](#)