



Product description

Order number	86.1171.001
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl, injected graduation up to 3 ml, sterile, 1 piece(s)/peel-pack

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3,5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product in accordance with MDD or IVDD IVD product
CE certified CE
Purity standard sterile
Sterilisation Gas sterilisation
Lot or serial number yes

Packaging

Minimum order qty. 840
Type of inner packaging peel-pack
Piece(s) / inner box 42
Piece(s) / case 840
Piece(s) / pallet 30240
Depth of case 578 mm
Width of case 368 mm
Height of case 195 mm
Case volume 0,0415 cbm
Weight of product 0,0015 kg
Weight of case 2,6 kg
EAN of inner box 4038917126267
EAN case 4038917020541

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-07-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT