

Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
Article - No:	12-521-0000	Revision: 1
Remarks		single use only

Specification

Product	Stirring spatula PS, Ø3x120mm, highly transparent
Dimensions	diameter: 3,00 mm (+/- 0,2 mm)
	length: 120,00 mm (+/- 0,2 mm)
Colour	highly transparent
Material Characteristics	polystyrene (PS)
Temperature range	-10°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 500
	pcs per carton: 10.000
	pcs per pallet: 1.080.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	