

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>09-021-5020</b>	<b>Revision: 6</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Petri dish PS, Ø90x16,2 mm, without vents, (non slippery / stackable), transparent, sterile R
<b>Dimensions</b>	diameter: 90,00 mm (+/- 0,2 mm)
	height: 16,20 mm (+/- 0,2 mm)
<b>Colour</b>	transparent
<b>Material Characteristics</b>	polystyrene (PS)
<b>Temperature range</b>	-20°C - 60°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
<b>Purity Grade</b>	do not contain any pyrogens
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Shelf - Life</b>	36 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 20
	pcs per carton: 480
	pcs per pallet: 9.600
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

