

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	09-231-0000	Revision: 3
Remarks		single use only

Specification

Product	Petri dish PS, Ø120x120x16 mm, with 4 vents, transparent
Dimensions	length: 120,00 mm (+/- 0,2 mm)
	width: 120,00 mm (+/- 0,2 mm)
	height: 16,00 mm (+/- 0,2 mm)
Colour	transparent
Material Characteristics	polystyrene (PS)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 14
	pcs per carton: 252
	pcs per pallet: 3.024
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	