

## P-9/500, platform

### DESCRIPTION

Platform with clamps for flasks.

### SPECIFICATIONS

Working area	360 × 400 mm
Flask volume	500 ml
Number of clamps	9



### CAT. NUMBER

BS-010135-AK

Platform P-9/500

### ACCESSORIES



ES-20/60  
Without platform  
Orbital Shaker-Incubator

Orbital Shaker-Incubator ES-20/60 for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



PSU-20i  
Without platform  
Multi-functional Orbital Shaker

Shaker PSU-20i provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)