

Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
Article - No:	09-385-0000	Revision: 1
Remarks		single use only

Specification

Product	Security closure PA, for sterilisation double bag 09-382-0000, natural, autoclavable
Colour	natural
Material Characteristics	polyamide (PA)
Autoclavability	autocl.
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Packaging	pcs per bag: 1.000
	pcs per carton: 1.000
	pcs per pallet: 117.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	