

Technical Datasheet

EN ISO 13485 certified	
Article - No:	12-521-9101
Remarks	Revision: 4 single use only

Specification

Product	Aspirating pipette PS, 2ml, single peel-pack, transparent, np pcr ready, sterile R
Dimensions	length: 275,00 mm (+/- 0,2 mm)
	diameter: 6,22 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	2ml
Material Characteristics	polystyrene (PS)
Temperature range	room temperature up to max. +37°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 50
	pcs per carton: 500
	pcs per pallet: 50.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

