

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-272-2100	Revision: 1
Remarks		single use only

Specification

Product	Mailing container PE, Ø29x85 mm, for 2 microscope slides, natural
Dimensions	inner diameter opening: 29,00 mm (+/- 0,2 mm)
	height: 85,00 mm (+/- 0,2 mm)
Colour	natural
Material Characteristics	polyethylene (PE)
Temperature range	-50°C – 70°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 500
	pcs per carton: 500
	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

