

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 11-652-0000 | Revision: 2 |
| Remarks | | single use only |

Specification

| | |
|--|--|
| Product | Screw cap container PP, 100ml, Ø48x80 mm, with enclosed screw cap PE, container transparent, cap natural |
| Dimensions | diameter: 48,00 mm (+/- 0,2 mm) |
| | height: 80,00 mm (+/- 0,2 mm) |
| Colour | container transparent, cap natural |
| Nominal Capacity | 100ml |
| Material Characteristics | container: polypropylene (PP) |
| | cap: polyethylene (PE) |
| Temperature range | -10°C - 40°C |
| Product Characteristics | biologically inert |
| | physiologically harmless |
| | no heavy metals according to EC/94/62 |
| | flexible and unbreakable |
| Shelf - Life | not limited at controlled storage conditions |
| Packaging | pcs per bag: 250 |
| | pcs per carton: 250 |
| | pcs per pallet: 6.000 |
| Storage | dry condition: 40%-60% rel. humidity |
| | bags / cases orig. closed |
| | storage temp. +10°C to max. +30°C |
| | pal. not stackable |
| Approved by nerbe plus quality department | |

00-000-0001/06.2020

