

# Technical Datasheet

Date: 13.10.23

Page 1 of 1

|                        |             |                 |
|------------------------|-------------|-----------------|
| EN ISO 13485 certified |             |                 |
| Article - No:          | 08-022-6004 | Revision: 2     |
| Remarks                |             | single use only |

## Specification

|                                           |                                                                                                                                                                   |
|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product                                   | Universal specimen container with spoon PP, 15ml, Ø20x76mm, with fitted screw cap PE, skirted, container transparent, cap brown, grad., sterile R max. RCF 2.500g |
| Dimensions                                | diameter: 20,00 mm (+/- 0,2 mm)                                                                                                                                   |
|                                           | height (without screw cap): 76,00 mm (+/- 0,2 mm)                                                                                                                 |
| Colour                                    | container transparent, cap brown                                                                                                                                  |
| Graduation                                | grad.                                                                                                                                                             |
| Nominal Capacity                          | 15ml                                                                                                                                                              |
| Material Characteristics                  | tube: polypropylene (PP)                                                                                                                                          |
|                                           | cap: polyethylene (PE)                                                                                                                                            |
| Temperature range                         | -20°C - 65°C                                                                                                                                                      |
| Centrifugation                            | max. RCF 2.500g                                                                                                                                                   |
| Product Characteristics                   | biologically inert                                                                                                                                                |
|                                           | physiologically harmless                                                                                                                                          |
|                                           | no heavy metals according to EC/94/62                                                                                                                             |
| Purity Grade                              | do not contain any pyrogens                                                                                                                                       |
| Sterilisation Process                     | radiation sterilization SAL 10-6                                                                                                                                  |
| Packaging                                 | pcs per bag: 100                                                                                                                                                  |
|                                           | pcs per carton: 1.000                                                                                                                                             |
|                                           | pcs per pallet: 24.000                                                                                                                                            |
| Storage                                   | dry condition: 40%-60% rel. humidity                                                                                                                              |
|                                           | bags / cases orig. closed                                                                                                                                         |
|                                           | storage temp. +10°C to max. +30°C                                                                                                                                 |
|                                           | pal. not stackable                                                                                                                                                |
| Approved by nerbe plus quality department |                                                                                                                                                                   |