

Technical Datasheet

Date: 21.11.25

Page 1 of 1

EN ISO 13485 certified		
Article - No:	09-613-5043	Revision: 5
Remarks		single use only

Specification

Product	Blender bags PE, 400ml, 190x300 mm, lateral filter, 50 pcs/bag, transparent
Dimensions	width: 190,00 mm (+/- 0,2 mm)
	length: 300,00 mm (+/- 0,2 mm)
	thickness: 70 µm
Colour	transparent
Nominal Capacity	400ml
Material Characteristics	bag: polyethylene (PE)
	filter: polyester fleece
Temperature range	-40°C - 80°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	filter porosity < 250 microns
	free of detectable phthalate
Purity Grade	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months
Packaging	pcs per bag: 50
	pcs per carton: 500
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	